Rosetta[®] SM

Instructions for Use









Rosetta[®] SM

Mechanical properties

Mechanical tests on lithium disilicate (LS2) glass-ceramics after final heat-treatment show that Rosetta SM crowns can achieve a higher strength (440 MPa) which is up to that of existing commercial product.


Machining performance

Machinable lithium metasilicate (LS) glass-ceramics allow for excellent edge stability. Rosetta SM chairside blocks with LS crystals can achieve a lower biaxial strength, which enables easier wet grinding compared to existing commercial product.

Table of Contents

1	Indications / Preps Guide	;
2	Select the Blocks	
3	Imaging	
4	Milling	
5	Preparing for Crystallization	
6	Characterizing	
7	Preparing for Cemention	

1 Indications / Preps Guide



2 Indications



Inlays



Onlays



Veneers



Anterior Single Crowns




Posterior Single Crowns



3-Unit Bridge
*up to the second Premolar

3 Imaging*



4 Milling*



5 Preparing for Crystallization





Finish surface with griding instrument.*

Use the fix material* before heat-treatment for crystallization.

6 Characterizing




Crystallization Schedule


Entry Temperature	Heating Rate	Final Temperature	Holding Time	Lowering Table	Vacuum On	Vacuum Off
400 °C	60 °C / min	840 °C	10:00 min.	700 °C	550 °C	840 °C

① Note

There may be a slight difference between the displayed temperature and the actual temperature of each furnace. Before you sinter Rosetta [®] SM blocks, please verify that the above recommended schedule is suitable for the furnace being used. Otherwise, try to find the optimized crystallization temperature through the following process.

- If there are some changes in the shape of restoration during crystallization heat treatment, please reduce the final temperature by 10~15° C.
- * Please follow our recommended crystallization schedule for best results.





O Note

Layer thickness

Dimension in mm

А	1.0	1.5	2.0	2.5	3.0
В	0.5	0.8	1.1	1.3	1.6
С	0.5	0.7	0.9	1.2	1.4


- A: Overall thickness
- B : Framework thickness
- C: Veneering material* thickness

B > C, Framework is to be thicker than veneering material*

7 Preparing for Cementation



Do not blast restoration.



Etch for 20 sec with 5% hydrofluoric acid.*

^{*} Respect all information given in the manufacturer's usage regulations.