

EN



Lithium Disilicate-Based Press on Zirconia

# Amber<sup>®</sup> LiSi-POZ

User's Manual



[www.hassbio.com](http://www.hassbio.com)



RX Only




Human-Aid  
System Supplier


Amber<sup>®</sup> LiSi-POZ User's Manual

Table of Contents

1 Introduction	3
2 Preparation guide	4
3 Zirconia framework	5
4 Wax-up	9
5 Sprueing	10
6 Ingot Selection	11
7 Investing	12
8 Burn-Out	13
9 Heat-Pressing	14
10 Divesting	15
11 Characterizing	16
12 Staining & Glazing	17
13 Completion	18
14 Indications / Contra-Indications	19
15 Product Line-up	20



## 2. Preparation Guide




... The minimum thickness of zirconia framework should be more than 0.6mm

### TIP!

- ! Make the prep tooth surface in the most rounded shape possible.(Deep chamfer margin, rounded shoulder margin).
- ! Maintain the most even margin thickness possible.


## 3. Zirconia framework

- ... Once the restoration shape is finished in CAD software, complete zirconia framework design with the cut-back technique.
- ... In selecting the shade for zirconia framework, please select one step brighter one than the shade you planned for the final restoration.





### TIP!

- ! Please keep the minimum thickness of cut-back to 0.6mm



- ... Trim the surface of zirconia framework so that edge areas would be a rounded shape instead of sharply angled.
- ... Keep the margins clean, not leaving the surface rough and get the space for inserting Amber LiSi-POZ.




- ... Zirconia framework should be treated by sandblasting, which entailed spraying 50~80µm of alumina( $Al_2O_3$ ) powder at a pressure of 2 bar and sintered at 1050°C for 15 minutes to stabilize zirconia framework.

**TIP!**

! Thermal stabilizing ※ The schedule table below is for HASS Zirtooth. Please refer to the firing schedule for each brand.

Starting Temp.	Up / Dry	Heating Rate	Max Temp.	Holding Time	Vacuum
500°C	1 min	65°C/min	1,050°C	15 min	NO

! Under cut design of framework creates mechanical adhesion which enhances the bond strength.



Incisial edge



Occlusal edge




Grooved surface

## 4. Wax-up

### TIP!


! In the anterior, pay attention to the occlusal surface during design.



Good




Acceptable





Not recommended

! For the stability of press restorations, consider 1~1.5mm width in designing.



! For decreasing tensile strength, make sure the ledge height on the lingual side is more than 3mm.




... After firing, measure the overall weight of zirconia framework. (It is important to decide the size of ingots)

... Complete the final shape of restorations. Remember to use combustible wax when doing a burn-out process.



Anterior Crown




Posterior Crown

### TIP!

! Form shapes while ensuring the wax thickness is not less than 0.6mm

## 5. Sprueing

Attach the sprues in the direction of flow for ceramic so that ingot can flow smoother during pressing.



Connect the object and investment ring base at an  $\angle 45\sim 60^\circ$  angle, at a length of 3~8mm, using 3~3.5 mm of spruing wax.



Keep a distance of at least 5 mm between the wax-up objects and silicone ring.

It is recommended to attach sprueing wax to each crown and it aids gas ventilation if air vent is attached in the thick part.


When finishing sprueing work, measure the total weight and subtract the weight of zirconia framework to decide the proper ingot size.



Ingot	Wax Weight	Invest. Ring
R10 1 ea(3 g)	~ 0.7 g	100 g
R15 1 ea(4.5 g)	0.7 ~ 1.2 g	200 g
R20 1 ea(6 g)	1.2 ~ 1.4 g	200 g

## 6. Ingot Selection

Color Chart



Conversion Chart

		Amber® LiSi-POZ				
		LT 0 / HT 0 (W4)	LT 1 / HT 1 (A1)	LT 2 / HT 2 (A2)	LT 3 / HT 3 (A3)	LT 4 / HT 4 (A3.5)
Vita Classic Shade	BL	BL1 / BL2				
	A	A1	A1 / A2	A2 / A3	A3.5	A4
	B	B1	B2	B2 / B3	B4	B4
	C		C1	C2 / C3	C3	C4
	D			D2 / D3	D3	D3

**TIP!**


Please choose one step brighter shade than the one you actually plan for the final restoration. (This prevents restoration from turning greyish during staining.)

Please choose the ingot which is most closely matched with the adjacent teeth.



## 7. Investing

- After mixing powder and liquid by hand for 20 seconds, mix it again with vacuum mixer. If it has hardened in the pressurizer after investing, strength and surface roughness are enhanced during pressing.



### TIP!

- For details, please refer to the IFU from the investment material manufacturer.

Phosphate-based investment material for ceramic press

# Amber<sup>®</sup> Vest





Packaging: KIT POWDER + EXPANSOR

**Amber<sup>®</sup> Vest POWDER**  
5kg (50X100g)

+

**Amber<sup>®</sup> Vest EXPANSOR-B**  
LIQUID (1,000ml)



Comparison of Reaction Layer Generation on Surface

## 8. Burn-Out




- Remove the silicone ring only after the investment is completely set.
- Trim the upper side flat and place the investment ring in the preheating furnace.
- The lower side of the investment should face down. Pay attention to ensure good drainage of the melted wax.

### TIP!


- Burn-out temperature and time should be according to the manufacturer's guidelines.

ex) Phosphate-based investment material for ceramic press **The highest temperature : 850°C**


# Amber<sup>®</sup> Vest



## 9. Heat-Pressing



... Make sure to put the ingot and plunger into the ring only at room temperature. At this time, printed side of the ingot should face up. Check if the ring bottom is placed flat.



... Proceed to pressing the ingot at the appropriate temperature.

### ... Pressing Schedules

	Translucency	Size	Shade	Investment Ring	Starting Temp.	Heating Rate	Max Temp.	Holding Time	Vacuum On	Vacuum Off
*Horizon	HT	R10(3 g) / R15(4.5 g)	W4, A1, A2, A3, A4	Small (100 g)	700	45°C/min	915°C	15 min	700°C	915°C
	LT			Large (200 g)				30 min		
	HT	R20 (6 g)		Small (100 g)				15 min		
	LT			Large (200 g)				30 min		

\*Horizon is a registered trademark of Shenzhou.



	Translucency	Size	Shade	Investment Ring	Starting Temp.	Heating Rate	Max Temp.	Holding Time	Press duration	Press level
*Austromat Press-i-dent	HT	R10(3 g) / R15(4.5 g)	W4, A1, A2, A3, A3.5	Small (100 g)	700°C	45°C/min	930°C	20 min	Auto1	6
	LT			Large (200 g)				30 min		
	HT	R20 (6 g)		Small (100 g)				20 min		
	LT			Large (200 g)				30 min		

\*Austromat Press-i-dent is a registered trademark of DEKEMA.

### TIP!

- ! There may be a difference between the temperature indicated on the furnace and the actual temperature. If problems occur after pressing, find out the optimal pressing temperature with the following process.
  - Bubbles or discoloration on restoration surface : Decrease the max temperature by 5~10°C degrees and try again.
  - If pressing is not completed : Increase the max temperature by 5~10°C degrees and try additional 5 minutes of holding time.
- ! Please do not use two of R 10 ingots so that air trap problem would not happen. R20 ingot required.

## 10. Divesting



... First check the length of the plunger and cut the investment with a disk.

... Use Al<sub>2</sub>O<sub>3</sub> for sandblasting. 4 bar of pressure for general blasting and 2 bar for precise blasting is recommended. Be cautious and only work after the ring has fully cool down.

### TIP!

- ! When cutting sprues, keep getting disk wet with plenty of water so that you can be cautious about micro fracturing. Refer to the instructions for use of the corresponding investment materials. Just few amount of reaction layer remains on the result at the recommended temperature.




## 11. Characterizing



Trim sprue and additionally layered areas. Ensure the surface is clean by removing bubbles. At this time, it is also necessary to work while applying water.

## 12. Staining & Glazing



After contouring, sandblasting the spot with  $\text{Al}_2\text{O}_3$  where staining procedures would be done, using 1 bar or less pressure. Apply the stain in accordance with the target shade.

## 13. Completion

### Anterior



Courtesy of CDT. Won Pil Jang and Dr. Hee Kyong Lee, Seoul, Korea

### Posterior



Courtesy of CDT. Won Pil Jang and Dr. Hee Kyong Lee, Seoul, Korea

## 14. Indications / Contra-Indications

### Indication



Single crown  
Press over zirconia framework



Multi unit bridge(3-unit / 4-unit)  
Press over zirconia framework

Amber LiSi-POZ  
Zirconia framework



Press over implant abutment made of zirconia

### Contraindication

- Subgingival abutment




- Press over unsintered zirconia framework

- Bruxism

## 15. Product Line-up



Human-Aid System Supplier **HASS** *beLIVE*

Amber <sup>®</sup> LiSi-POZ		Dimensions (mm)	pcs / Pack
	R10	Ø12.7 × 10T	5 Ingots
	R15	Ø12.7 × 15T	3 Ingots
	R20	Ø12.7 × 20T	3 Ingots

### HASS Corporation

77-14, Gwahakdanji-ro, Gangneung-si, Gangwon-do, KOREA 25452  
 Tel: +82-70-7712-1300 / Fax: +82-33-644-1231  
 Customer Support : +82-2-2083-1367  
 E-mail : hasscorp@hassbio.com  
 Website : www.hassbio.com

Printed in KOREA © HASS Corporation. All rights reserved.

This material is designed for usage in dentistry. Follow instructions HASS is not liable for any loss caused by failure to comply with regulation or scope of indication. Users are responsible for testing products to verify the compatibility for any usage which are not written in the instructions. The explanations and data contained within do not carry any guarantees and/or obligations. All enclosed recommendations and restrictions apply when used with products from other manufacturers.