

EN

Phosphate-based investment material for ceramic press



User's Manual



Working instructions



1. Liquid/Cylinder Preparation

Concentration of mixing liquid

Indication	Liquid (EXPANSOR-B)	Distilled water
Inlays 🖷	13 ml (50%)	13 ml (50%)
Veneers @	15.6 ml (60%)	10.4 ml (40%)
Crowns $\overline{\mathbb{W}}$	18.2 ml (70%)	7.8 ml (30%)

IMPORTANT:

Amber® Vest must always be used with its special mixing liquid EXPANSOR-B, pure or diluted. The dilutions can be effected by adding distilled or demineralised water, shaking slightly and then applying to the bottle an identifying label of the concentration made.



DO NOT USE EXPANSOR LIQUID WHICH HAS BEEN FROZEN. (BECOME CRYSTALLISED)



2. Mixing



- · To obtain good results, it is essential to use the precise measurements.
- · For greater ease of use, Amber Vest is supplied in predosed sachets of 100g, and a dosing tube is supplied for the liquid.
- \cdot It is always recommended to use whole sachets and to measure the liquid with the original dosing tube.
- · In the case of using non-pre-dosed portions of powder, weigh the exact contents with a scale.

The powder liquid ratio 100g powder: 26 ml liquid



DO NOT USE BOWLS NEITHER CONTAINERS WHICH HAVE BEEN IN CONTACT WITH PLASTER.





- $\cdot \ \text{Prepare the liquid EXPANSOR-B according to the concentration required (as indicated above)}.$
- \cdot Add the powder to the liquid and mix vigorously by hand to obtain a totally moist and uniform mass.
- · Then pass to mechanical vacuum mixing for 60 seconds depending on the mixer.

3. Filling the cylinder



- · Before filling the cylinder, apply a degreasing agent or wax cleaner with a spray or a brush uniformly and dry for a few seconds. It may allow a better adherence and flows of the investment.
- · When the mixing time is completed, fill the cylinder, vibrating gently, for 10-15 seconds, then stop.



4. Working Time



6-7 minutes (20-22°C) from the moment that powder and liquid contact.

The working time can be reduced:

- · When the liquid is at a higher temperature
- · With ambient temperatures above 22°C
- · With the use of contaminated utensils

5. Burn out

Heating can be carried out in two different ways.

5-1. Rapid burn out

- 1. Before investment, preheat the ring furnace to the final temperature.
- 2. Proceed with the investment process.
- 3. Place the ring in the ring furnace 20-30 minutes after the investment, before it reaches the final temperature. Press at the final temperature (850°C)
 - * Caution: the furnace temperature may vary depending on the brand and environmental conditions.

5-2. Standard burn out

- 1. After allowing the cylinder to settle at room temperature, place it in the furnace and maintain a temperature of 250°C for 40 minutes.
 - * Caution: for a safer procedure, it is recommended to anchor the temperature at 500°C for approximately 30 minutes.
- 2. Gradually increase the temperature to the final target casting temperature at a rate 7°C/ min.



6. Final Result







* HS Code: 3407.00

* Classification: Type1, Class2



Manufacturer

Famadent, S.L.U. (PROTECHNO®)

Pol. Ind. Emporda Internacional, CI Garrotxa, 6 17469 Vilamalla, Girona. Spain
Tel.: (34) 972 52 61 69 / Fax: (34) 972 52 55 86
protechno@protechno.com / www.protechno.com

Distributor

HASS Corporation

77-14, Gwahakdanji-ro, Gangneung-si, Gangwon-do, KOREA 25452 Tel: +82-70-7712-1300 / Fax: +82-33-644-1231

Customer Support: +82-2-2083-1367 E-mail: hasscorp@hassbio.com Website: www.hassbio.com

HASS BIO America, Inc.

10400 Eaton place, Suite 220 | Fairfax, VA | 22030 Tel : 703-537-0333

E-mail: info@hassbioamrica.com Website: www.hassbioamerica.com This material is designed for usage in dentistry. Follow all instructions HASS is not liable for any loss caused by failure to comply with regulation or scope of indication. Users are responsible for testing products to verify the compatibility for any usage that is not written in the instructions. The explanations and data contained herein

do not carry any guarantees and/or obligations.

All enclosed recommendations and restrictions apply when used with products from other manufacturers.