

Capture the Truth

Natural Daylight Solution for Mobile
Dental Photography

COCO Lux[®]



www.hassbio.com

CE2195



RX Only

RoHS

HASS | Human-Aid
System Supplier

EASY & SIMPLE

Sick and tired of heavy, complicated DSLR package? COCO Lux[®] with your smartphone is the right solution for your dental photography.

EASY PHOTOGRAPHING


Just attach COCO Lux[®] onto the smartphone - and you are ready to go for pictures. COCO Lux[®] is compatible with smartphones that are 55 - 89mm wide and is able to change its position according to different camera location of each phone and it even fits phones with case.*

* Flip case may require removal of its cover to use this solution.



SIMPLE SHARING & COOPERATION

Diverse messenger apps on smartphone will enable you to experience instant and interactive communication simply after sharing the resulted pictures taken by COCO Lux[®].



TAKE

&



SHARE


TRUE & VERSATILE

In addition to its simplicity and ease, COCO Lux[®] also features the indispensable essence of photography - desirable light settings for photography and practical values for a wide variety of purposes.

GREAT ACCURACY & DEFINITION


COCO Lux[®] hires LED light source that is the closest to natural daylight: 5,500K color temperature (white) and over 90Ra color rendering index.* Polarizer is also bundled for the effective reflection management.

* Color rendering index (CRI) indicates how similar a specific light is to natural daylight in percentage.



HIGH VERSATILITY

As COCO Lux[®] supports both twin flash mode and ring flash mode, every user may make use of it for a number of different purposes such as taking shade, clinical case pictures, model work photograph, video, etc. Lighting part of COCO Lux[®] is rotatable in 90° range for twin flash mode both vertically and horizontally. Plus, either pair of the lighting part is movable separately. Use your imagination and creativity to discover unprecedented images and quality from COCO Lux[®].



Courtesy of CDT. Min Goo Cho, Seoul, Korea



Courtesy of DR. Sung-ki Kim, Busan, Korea

Natural Daylight Solution for Mobile
Dental Photography

COCO Lux[®]

HASS Corporation

77-14, Gwahakdanji-ro, Gangneung-si, Gangwon-do,
KOREA 25452
Tel: +82-70-7712-1300 / Fax: +82-33-644-1231
Customer Support : +82-2-2083-1367
E-mail : hasscorp@hassbio.com
Website : www.hassbio.com

HASS BIO America, Inc.

10400 Eaton place, Suite 220 I Fairfax, VA I 22030
Tel : 703-537-0333
E-mail : info@hassbioamerica.com
Website : www.hassbioamerica.com

This material is designed for usage in dentistry. Follow instructions HASS is not liable for any loss caused by failure to comply with regulation or scope of indication. Users are responsible for testing products to verify the compatibility for any usage which are not written in the instructions. The explanations and data contained within do not carry any guarantees and/or obligations. All enclosed recommendations and restrictions apply when used with products from other manufacturers.